



# Payload-Linker Services for Antibody Drug Conjugates

# ADC Payload & Linker Solutions

SK pharmteco has extensive experience in safely handling and manufacturing high-potency APIs (HPAPIs) and intermediates down to 10ng/m<sup>3</sup>. Capitalizing on this expertise, we provide custom development and CGMP manufacture of ADC payload-linkers from gram to kilogram scale.

## OUR FOCUS

### Delivering value to you



**Experience:** We support commercial programs with leading ADC developers to provide payload and linker components. We offer a seamless partnership with prompt communication and a dedicated program manager to keep you informed throughout your project.



**Expertise:** We have decades of experience in HPAPI manufacturing (to 10 ng/m<sup>3</sup>) with 50+ commercial batches produced. Additionally, as a leader in chromatography, we offer unparalleled support for the purification of payloads.



**Quality:** We adhere to the highest standards of quality, supported by a successful track record with the FDA, EMA, PMDA, and other agencies.



**Scalability:** Our scalable production capabilities support grams to kilograms, to accelerate your time-to-market.



**Proprietary Payload & Linker Fragments:** Our established licensing agreements provide broad access to proprietary fragments of payloads and linkers.



**Comprehensive Analytical Services:** In addition to our manufacturing capabilities, we offer comprehensive analytical testing services to ensure the highest standards of quality and compliance.

## OUR FACILITIES

### Flexibility Meets Capacity

#### GMP Analytical Laboratory

**Batch size:** < 0.2 kg

**Containment:** OEL < 10 ng/m<sup>3</sup>

**Equipment:**

- Up to 5 L vessels
- HPLC, TFF & freeze drying
- Drying capability qualified from mL-L
- Easymax (programmable reaction calorimeter) supports safety studies & band 5 materials
- High containment jet milling (2 inch)

**Operational**

#### Small-Scale GMP Facility

**Batch size:** 0.2 - 0.5 kg

**Containment:** OEL < 10 ng/m<sup>3</sup>

**Equipment:**

- 3 reactors (up to 30 L)
- Glovebox isolators (rigid or flexible)
- High pressure HPLC (up to 150 mm column)
- TFF filtration
- Stainless steel tray dryer (0.2 m<sup>2</sup>)
- Lyophilizer (ice capacity 12 kg)

**Coming online by Q1 2026**

#### GMP Manufacturing Plant

**Batch size:** > 0.5 kg

**Containment:** OEL < 10 ng/m<sup>3</sup>

**Equipment:**

- 3 reactors (60-90 L) & associated tanks (60 L EA)
- Glovebox isolators (rigid)
- Flash Chromatography (150 L column)
- Hastelloy 25 L / 0.3 m<sup>2</sup> filter dryer

**Operational**